



Laboratory Colour Gauge

With Automatic Turntable

CR100LG

CRI00LG

Applications Include:

Food Products

- Potato Crisps
- Corn Snacks
- Tortillas
- Sugar
- Peanuts
- Breakfast Cereals
- Pasta
- Biscuits
- Coffee

Chemicals/Minerals

- Detergents
- Ores
- Ceramic Powders



Laboratory Colour Gauge

White Light Colour Measurement

The **CRI00LG** uses the latest solid state white light LED to illuminate the product to be measured. A variety of colour space outputs are available such as the widely accepted CIE L* a* b* standard, or other common colour spaces. The **CRI00LG** Can be used for quality control in a wide variety of industrial processes , providing fast accurate results every time....

The **CRI00LG** Bench is rugged , simple to use and highly reliable on the factory floor or quality control lab.

CRI00LG Features:-

- Bright Touch Screen LCD above sample area - Less effected from sample debris
- Programmable Turntable automatically rotates
- Rugged Construction
- USB for data download to memory stick

Simple , Easy, Reproducible Results

- Fill sample dish to constant volume
- Place dish onto Turntable
- Dish automatically rotates
- Wait 10 seconds
- Read Colour value

Sample Preparation

Many products require no sample preparation. Sampling error is reduced because of the rotating dish that permits the **CRI00LG** to “see” multiple coloured surfaces and average hundreds of measurement readings.



Senware Limited

Unit 4, Adelaide House,
CorbyGate Business Park
Priors Haw Road
Corby, Northants
NN17 5JG

For More Information :-

Tel:- +44 (0)1536 408066
Fax :- +44 (0)1536 407813
Email :- sales@senware.co.uk
Website www.senware.co.uk